

# I N S T R U C T I O N S




## X-Fer Rail™ Folding Transfer Board

70404-0000

Before using the X-Fer Rail™ consult your therapist or physician regarding safe transferring techniques


### Unfolding X-Fer Rail™:

1. Unwind strap
2. Hold end of the X-Fer Rail™ and allow both halves of the rail to come together
3. Make sure both halves are fully engaged. X-Fer Rail™ is ready to use



### Folding X-Fer Rail™:

1. Pull halves apart to disengage connectors
2. Fold the halves of the X-Fer Rail™ until the ends come together
3. For compact storage wrap strap around the ends of the X-Fer Rail™



Weight Capacity: 250 Lbs. (113.6 Kg)

**Warning:** The elastic cord on the inside of the X-Fer Rail™ should be replaced every 2 years or whenever visible fraying occurs. Please contact Maddak Customer Service at [custservice@maddak.com](mailto:custservice@maddak.com) or 1-800-443-4926 to arrange replacement of the cord.

## Maddak

661 Route 23 South • Wayne, NJ 07470-6814 USA  
973-628-7600 • fax 973-305-0841 • email [custservice@maddak.com](mailto:custservice@maddak.com) • [www.maddak.com](http://www.maddak.com)