


Putting on Panties, Shorts, Pajama Bottoms, Slacks, etc.

While sitting, hold the two parts of SLIP-ON[®], one in each hand, and engage the hook end with the inside of the waistband of the apparel. Lower the apparel to the floor and pull it over the feet and up both legs until it can be reached by the hands and pulled into the desired position.

If the waistband has an opening, as in slacks, close it before beginning and reopen before pulling over the hips.



SLIP-ON[®]

CATALOG NO. 73839-0000

MAKES DRESSING EASY



Maddak Inc.

6 Industrial Road • Pequannock, NJ 07440

tel: 973-628-7600 • fax: 973-305-0841

email: custservice@maddak.com • www.maddak.com


SLIP-ON[®] helps you put on the following apparel without the aid of another person:

Stockings—Socks • Shorts
Panties—Undershorts
Pajama Bottoms • Slacks

Can also be used to remove shoes.

INSTRUCTIONS

Engage the gear-like teeth to make SLIP-ON® operable. As shown in Figure 1, hold SLIP-ON® in one hand and insert the hook ends inside the top opening of stocking. Hold SLIP-ON® in the left hand if you are putting the stocking on the left foot; in the right hand for the right foot. Squeeze to spread wide enough to firmly hold the stocking, with the heel hanging down.



Geer Teeth




Fig. 1

As shown in Figure 2, continue squeezing to spread the stocking top wide enough to slip over the toes and pull it back toward the heel. Inching the foot along the floor while pulling is helpful.



Fig. 2

As shown in Figure 3, keep the hook ends engaged with the stocking and separate the two parts of SLIP-ON®, holding one in each hand. Pull the stocking over the heel and up the leg to the point where it can be reached by hand and adjusted to the desired position.

Caution: When using with sheer stockings engage hook ends with reinforced top band to avoid rupturing the sheer knit.




Fig. 3

Removing the Shoe

Engage the notch end of one of the parts of SLIP-ON® with the counter of the shoe and push down—at the same time raise the heel slightly and the shoe comes off easily. Little pressure is required, especially with slip-on shoes.



Fig. 4