

I N S T R U C T I O N S



Basic Elevated Toilet Seat

Catalog Number 72579-0000 Basic Elevated Toilet Seat
Catalog Number 72579-0001/0011 Basic Elevated Toilet Seat
with Closed Front Adapter

1. Lift existing toilet seat and cover
2. Place the Basic Elevated Toilet Seat on top of the toilet. Make sure the flange of the elevated toilet seat fits within the rim of the toilet.
3. For closed front models, simply slip the adapter over the front of the toilet



Maddak

6 INDUSTRIAL ROAD • PEQUANNOCK, NJ 07440-1993 USA

973-628-7600 • fax 973-305-0841 • email custservice@maddak.com • www.maddak.com