

I N S T R U C T I O N S



Pistol-Grip Remote Toenail Clipper

Catalog Number 74139-0000

U.S. Patent #4,893,406

The Remote Toenail Clipper gives you an extra 18½" (47 cm) of reach that makes it possible to clip toenails without bending.

Before using, observe the operation of the clipper by squeezing the trigger and noting the motion of the clipping head.

To clip toenails, slip the clipping blades over the toenail and squeeze the trigger. Be careful not to pinch any part of the toe or to clip too deeply.

Release trigger and continue on.

Keep out of reach of children.



Maddak Inc.

6 INDUSTRIAL ROAD • PEQUANNOCK, NJ 07440-1993 USA

973-628-7600 • fax 973-305-0841 • email custservice@maddak.com • www.maddak.com