

I N S T R U C T I O N S



Press-On Nail Clipper

74140-0000

Cutting Finger Nails (excluding the thumb)

1. Place finger nail in between the two blades of the nail clipper
2. Press down on the platform with the palm of your hand, using enough pressure to cut through the nail



Cutting Thumb Nail

1. Place thumb nail in between the blades of the nail clipper
2. Create a fist
3. Push down on the platform with your fist, using enough pressure to cut through the nail



Maddak

6 INDUSTRIAL ROAD • PEQUANNOCK, NJ 07440-1993 USA
973-628-7600 • fax 973-305-0841 • email custservice@maddak.com • www.maddak.com