



Elevated Toilet Seat Dimensions Chart

ELEVATED PUSH-UP TOILET SEAT, 725600050



BOTTOM VIEW



SIDE VIEW

Item #	Description	A		B		C		D		E		F		G		H	
		Inches	CM	Inches	CM	Inches	CM	Inches	CM	Inches	CM	Inches	CM	Inches	CM	Inches	CM
725600050	Elevated Push-Up Toilet Seat, 3" (7.62 cm)	12	30.5	11	27.9	14	35.6	9	22.9	17 ³ / ₄	45.1	24 ¹ / ₄	61.6	1 ¹ / ₂	1.3	3	7.6



Line drawing is a general representation and may not actually portray the shape of this specific toilet seat. Please refer to photo and dimensions chart for specific information on this model to determine compatibility of elevated toilet seat with your existing toilet seat bowl.