



# Elevated Toilet Seat Dimensions Chart

TALL-ETTE ECONOMY W/BOLT-DOWN, 725812001



**BOTTOM VIEW**



**SIDE VIEW**

Item #	Description	A		B		C		D		E		F		G		H	
		Inches	CM	Inches	CM	Inches	CM	Inches	CM	Inches	CM	Inches	CM	Inches	CM	Inches	CM
725812001	Tall-Ette Economy W/Bolt-Down	10	25.4	7 1/2	19.1	11 1/2	29.2	9 1/2	24.1	14 1/4	36.2	16	40.6	1 1/2	3.8	4	10.2



Line drawing is a general representation and may not actually portray the shape of this specific toilet seat. Please refer to photo and dimensions chart for specific information on this model to determine compatibility of elevated toilet seat with your existing toilet seat bowl.