

I N S T R U C T I O N S



Reacher Clip

76972-1000

Keeps reachers handy and firmly attached to wheelchair and walker frames

The larger side of the double clip snaps onto the wheelchair, walker frame or bed rail and the smaller side clips onto a reacher. The Reacher Clip fits most reachers and can be snapped onto tubing at any angle and rotated to the correct position.



Maddak

6 INDUSTRIAL ROAD • PEQUANNOCK, NJ 07440-1993 USA

973-628-7600 • fax 973-305-0841 • email custservice@maddak.com • www.maddak.com