



One-Handed Canister Set

Easy to open, no twisting or turning required

- Ideal for children, seniors or anyone with limited hand strength
- Canisters open and close with a press of the finger or palm of the hand, requiring minimal force and strength
- Top rack dishwasher safe



Made in Denmark

Three sizes:

- Small Canister:
25 oz. (750 ml)
- Medium Canister:
45 oz. (1350 ml)
- Large Canister:
67 oz. (2000 ml)

Item 753600000



SP Ableware
A Division of SP Industries
www.spableware.com

Instructions:

1. To open simply press the small dot on the rim of the lid with your finger or the palm of your hand.
2. Lift to remove lid; no twisting or turning is required.
3. To close simply press down on the middle of the lid with your finger or the palm of your hand.

Note: The one-handed canisters are meant for storage purposes only – do not use in the microwave or freezer.

1



2



3



 Maddak
SP ABLEWARE
www.spableware.com