



Elevated Toilet Seat Dimensions Chart

HINGED TOILET SEAT, STD, 725711000



BOTTOM VIEW



SIDE VIEW

Item #	Description	A		B		C		D		E		F		G		H	
		Inches	CM	Inches	CM	Inches	CM	Inches	CM	Inches	CM	Inches	CM	Inches	CM	Inches	CM
725711000	Hinged Toilet Seat, Std	11¼	28.6	9½	24.1	11½	29.2	10¼	26.0	17½	44.5	14¼	36.2	½	1.3	3	7.6



Line drawing is a general representation and may not actually portray the shape of this specific toilet seat. Please refer to photo and dimensions chart for specific information on this model to determine compatibility of elevated toilet seat with your existing toilet seat bowl.