



# Elevated Toilet Seat Dimensions Chart

BOLT-ON TOILET SEAT W/ARMS, STD, 725753210



**BOTTOM VIEW**



**SIDE VIEW**

Item #	Description	A		B		C		D		E		F		G		H	
		Inches	CM	Inches	CM	Inches	CM	Inches	CM	Inches	CM	Inches	CM	Inches	CM	Inches	CM
725753210	Bolt-On Toilet Seat W/Arms, Std	11¼	28.6	9½	24.1	14⅞	37.8	6¼	15.9	17¼	43.8	21½	54.6	¼	0.6	3½	8.9

Line drawing is a general representation and may not actually portray the shape of this specific toilet seat. Please refer to photo and dimensions chart for specific information on this model to determine compatibility of elevated toilet seat with your existing toilet seat bowl.